



BEST ELECTRO
SURGICAL



Collection of Electro Surgical Instruments



Adalat Garh, Tehsil & District
Sialkot 51310 Pakistan +92-
329-6320000



info@bestelectrosurgical.com



www.bestelectrosurgical.com





BEST ELECTRO
SURGICAL



Adalat Garh, Tehsil & District
Sialkot 51310 Pakistan +92-
329-6320000



info@bestelectrosurgical.com
www.bestelectrosurgical.com

Introduction

Best Electro Surgical family has manufactured high quality surgical, dental, manicure, and veterinary instruments to the world market for over a century. More than 30 years ago, **Best Electro Surgical** established its own sales and marketing facilities in the Pakistan. Through constant innovation, we continue to provide the highest quality instruments and unparalleled service to meet and exceed all of our customer's needs.

Our customer base spans across the Pakistan, as well as globally. Although we have an extensive

customer base, we continue to treat each customer with individual attention. We are a smaller firm

that deals direct with manufacturer and therefore can offer competitive pricing and take pride in

giving a more personalized service. We view ourselves as an extension of your department and will

go to great lengths to find the instruments you require. We will provide competitive pricing on any

brand or manufacturer. We will do our best to work with your budget, and inventory any items upon

request.

Among surgical, dental, manicure, and veterinary instruments, sell a variety of surgical, dental,

manicure, and veterinary instruments here at **Best Electro Surgical**. We have handpicked some

of the industry's top products designed to relieve pain, increase range of motion, and get people

back on their feet as quick as possible. These products are used by physical therapists,

occupational therapists, athletic trainers, chiropractors, and the average person daily all over the

world.

We primarily distribute German manufactured, German stainless steel surgical, dental, manicure,

and veterinary instruments. All instruments have a LIFETIME GUARANTEE on material workmanship.

Thanks and regards,

Process of Instruments

Steel the quality of steel is very important in the production of quality instruments. The best you select the best you produce. In order to ensure the finest of quality steel, we have imported or best quality local steel as required by valued customers.

Forging is a process of cutting shaped raw instruments with the help of dyes and hammers. Here, the skilled hands give life to the metal and give it the rudimentary shape that will go through series of complicated processes.

Forging Inspection at this stage, forged tools are checked to ensure no piece contains any cracks or is reshaped beyond acceptable limits. Checking the quality of raw instruments for: -
a. Shape fault b. Broken

Machining in this process forging is cut or shaped by machines i.e. cutting sharp edges and raw material in the shape. Here instruments are given their shape through grinding and hammers. For this purpose, extremely skilled labor is required and this process is carried out by instrument specialists only.

Milling drilling of holes in the instruments where required for screws. Different kinds of screwing machines are used and instruments are made pairs here. Instruments with even the slightest difference in whole positions will suffer otherwise.

Inspection at this stage, all the instruments is inspected not only for machining and milling, but for every process they've been through. It's the stage where only best instruments are selected and instruments with tiny holes and manufacturing faults are again sorted out.

Filling at this stage the steel is not too hard so filling is used to shape the instruments in the exact shape they are required to be. The workers match the shape exactly with the samples provided at this stage. So the instruments get most of their shapes at this stage.

Filling Inspection after filling the inspection is made and measurements and shapes are checked according to samples.

Binding and Temper at this step the instruments are bound in groups of 12 pieces. Then they put them in temper machine, which is filled with a particular material. Here instruments are tempered and its ensured that every instrument get equal heat and treatment so a consistent quality would achieve.

Snuffing Fitting for scissors and forceps they are in two pieces so they are joined in this process with screw. Forceps (Tweezers) are joined before and their joint is leveled with a machine by burning the raw material on the joint. At the same stage, instruments are filled to eliminate welding suspensions.

Plant or Heat Treatment in this process the instruments are dipped in a material by hanging for five minutes to harden the material.

Polishing instruments is polished at this stage. Till this point at manufacturing stage, all the necessary shaping is done. Polishing makes every millimeter of instrument clear and any instrument with hidden defects are revealed could be sorted out.

Process of Instruments

Best Electro Surgical Inspection 2% of each lot of instruments is checked by executives and then goes the further processing. In inspection they see each and every part and defective parts are marked with red marker and if rejected send to polishing again.

Ultrasonic cleaning in this process the instruments are dipped in two materials for cleaning the polish raw thing on it. One min. for first material and then for half min. for other material. Ultrasonic cleaning cleans any polishing material stuck and leaves the instrument crystal clear.

Lubrication all the instruments involving screws are lubricated at this step. Lubricating the instruments moveable portions facilitates the instrument movement.

Tuck checking In this process workers check the cuts and all other things just by viewing the instruments individually.

Etching putting each instrument on stamp pad and burning process to put stamp on each instrument. This stamping procedure is done through electronic etching machines that make the stamp reliable and even sustainable under high boiling temperatures when instruments are sterilized.

Cleaning and Checking instruments are cleaned and checked again one by one. Again it's the checking for operation, measurements, polish and usability.

Best Electro Surgical Inspection 4 % of the instruments is inspection by SSE Directors. Samples for inspection are taken randomly and are observed to every bit of detail. Even the slightest difference from acceptable standards could result in rejection of whole lot.

Packing in boxes each contains 100 pieces. 10 Boxes in an each Carton. Here instruments are packed with respect to sizes and types. Instruments boxes are properly labeled and stored in proper storing environment.

Labeling is put in this process.

Storage by lot and order number

Table of Contents

01 European Bipolar Forceps	07
02 European Irrigation Bipolar Forceps	08
03 European Non-Stick Bipolar Forceps	08
04 Us Bipolar Forceps	11
05 Monopolar Forceps	17
06 Us Irrigation Bipolar Forceps	18
07 US Non-Stick Bipolar Forceps	19
08 Hand-Switch Bipolar Foreceps	21
09 Hand-Switch Monopolar Foreceps	22
10 Diathermy Instruments	23
11 Disposable Bipolar Forceps	25
12 OB/ Gyn Leep Instruments	27
13 OB/ Gyn Laser Surgery Instruments	30
14 Laparoscope & Arthsopice Electrode	32
15 Electro Surgical Pencil & Electrode	34
16 Electro Surgical Cable	39
17 Bipolar artery Sealer	42



BES-100
McPherson



BES-101
McPherson Ang.



BES-103
Jeweller #5 Str




BES-104
Jeweller #5 Ang



BES-105
Jeweller #7



BES-106
Adson



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP-** Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-100	McPherson	0.5 0.5	10.2
BES-101	McPherson Angled	0.5,	10.2
BES-102	Jeweller #5	0.7	11.5
BES-103	Jeweller #5 Angled	0.5,	11.5
BES-104	Jeweller #7	0.7	11.5
BES-105	Adson	0.5,	12.0
BES-106	Adson	0.7 1.0	15.0

BES-106
1.0
BES-
108
BES-

BES-114
Straight



BES-114-N
Straight



BES-116
Straight




BES-118
Straight



BES-120
Straight



BES-122
Straight



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-114	Straight	0.5, 1.0, 2.0	15.0
BES-114-N	Straight	0.5, 1.0, 2.0	15.0
B E S -	Straight	0.5, 1.0, 2.0	18.0
1 1 6	Straight	0.5, 1.0, 2.0	20.0
B E S -	Straight	0.5, 1.0, 2.0	22.0
1 1 8	Straight	0.5, 1.0, 2.0	24.0

B E S -

1 2 0

B E S -

122

BES-124
Straight Angled



BES-126
Straight Angled




BES-128
Straight Angled



BES-130
Straight Angled



BES-132
Straight Angled



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Straight	Tip Size (mm)	Length Size (cm)
BES-	Angled		0.5, 1.0, 2.0	15.0
124	Straight Angled		0.5, 1.0, 2.0	18.0
BES-	Straight Angled		0.5, 1.0, 2.0	20.0
126	Straight Angled		0.5, 1.0, 2.0	22.0
BES-	Straight Angled		0.5, 1.0, 2.0	24.0

128

BES-

130

BES-

132

BES-134
Bayonet



BES-134-N
Bayonet



BES-136
Bayonet



BES-138
Bayonet




BES-140
Bayonet



BES-142
Bayonet



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-134	Bayonet	0.5, 1.0, 2.0	16.5
BES-134-N	Bayonet	0.5, 1.0, 2.0	16.5
B E S -	Bayonet	0.5, 1.0, 2.0	18.0
1 3 6	Bayonet	0.5, 1.0, 2.0	20.0
B E S -	Bayonet	0.5, 1.0, 2.0	22.0
1 3 8	Bayonet	0.5, 1.0, 2.0	24.0

B E S -

1 4 0

B E S -

142

BES-144
Bayonet Angled Up



BES-146
Bayonet Angled Up




BES-148
Bayonet Angled Up



BES-150
Bayonet Angled Up



BES-152
Bayonet Angled Up



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP-** Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name Bayonet	Tip Size (mm)	Length Size (cm)
BES-	Angled Up	0.5, 1.0, 2.0	16.5
144	Bayonet Angled Up	0.5, 1.0, 2.0	18.0
BES-	Bayonet Angled Up	0.5, 1.0, 2.0	20.0
146	Bayonet Angled Up	0.5, 1.0, 2.0	22.0
BES-	Bayonet Angled Up	0.5, 1.0, 2.0	24.0

148

BES-

150

BES-

152

BES-154
Bayonet Angled Down



BES-156
Bayonet Angled Down




BES-158
Bayonet Angled Down



BES-160
Bayonet Angled Down



BES-162
Bayonet Angled Down



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-	Bayonet Angled Down	0.5, 1.0, 2.0	16.5
154	Bayonet Angled Down	0.5, 1.0, 2.0	18.0
BES-	Bayonet Angled Down	0.5, 1.0, 2.0	20.0
156	Bayonet Angled Down	0.5, 1.0, 2.0	22.0
BES-	Bayonet Angled Down	0.5, 1.0, 2.0	24.0

158

BES-

160

BES-


162



BES-170-IR-N
Bayonet Irrigation



BES-170-IR
Bayonet Irrigation




BES-168-IR
Bayonet Irrigation



BES-164-IR
Straight Irrigation

BES-166-IR
Straight Irrigation



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-164-IR	Straight Irrigation	0.5, 1.0, 2.0	18.0
BES-166-IR	Straight Irrigation	0.5, 1.0, 2.0	20.0
BES-168-IR	Bayonet Irrigation	0.5, 1.0, 2.0	18.0
BES-170-IR	Bayonet Irrigation	0.5, 1.0, 2.0	20.0
BES-170-IR-N	Bayonet Irrigation	0.5, 1.0, 2.0	20.0

Note: More Bipolar Irrigation Forceps are available use IR extension to any regular bipolar forceps.

BES-172-NS
McPherson



BES-174-NS
McPherson Angled




BES-176-NS
Jeweler # 5



BES-178-NS
Jeweler # 5 Angled



BES-180-NS
Adson



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-172-NS	McPherson	0.5,	10.2
NS	McPherson Angled	0.5,	10.2
NS	Jeweller #5	0.5,	11.5
BES-178-NS	Jeweller #5 Angled	0.7	11.5
BES-180-NS	Adson	0.5, 1.0	12.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.

BES-181-NS-N
Straight



BES-182-NS-N
Straight




BES-184-NS
Straight



BES-186-NS
Straight



BES-188-NS
Straight



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP-** Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-181-NS-N	Straight	0.5, 1.0, 2.0	15.0
BES-182-NS	Straight	0.5, 1.0, 2.0	15.0
BES-184-NS	Straight	0.5, 1.0, 2.0	18.0
BES-186-NS	Straight	0.5, 1.0, 2.0	20.0
BES-188-NS	Straight	0.5, 1.0, 2.0	22.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.

BES-190-NS
Bayonet



BES-190-NS-N
Bayonet




BES-192-NS
Bayonet



BES-194-NS
Bayonet



BES-196-NS
Bayonet




Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-190-NS	Bayonet	0.5, 1.0, 2.0	16.5
BES-190-NS-N	Bayonet	0.5, 1.0, 2.0	16.5
BES-192-NS	Bayonet	0.5, 1.0, 2.0	18.0
BES-194-NS	Bayonet	0.5, 1.0, 2.0	20.0
BES-196-NS	Bayonet	0.5, 1.0, 2.0	22.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.



BES-230
McPherson



BES-230-N
McPherson



BES-232
McPherson Ang.




BES-234
Coaptation



BES-236
Coaptation



BES-238
Coaptation



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-230	McPherson	0.5	10.2
BES-230-N	McPherson	0.5	10.2
BES-232	McPherson Angled	0.5	10.2
BES-234	Coaptation Nedler	0.5	10.2
BES-236	Coaptation	0.5	10.8
BES-238	Jeweller #5	Needle, 0.5	11.5



BES-240
jeweler # 4



BES-242
Tenzel




BES-246
iris Straight



BES-244
jeweler # 7





BES-248
iris Angled



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Jeweller	Tip Size (mm)	Length Size (cm)
BES-240	#4		Needle, 0.5	11.5
BES-242	Tenzel		Needle, 0.5	10.2
BES-244	Jeweller #7		Needle, 0.5	11.5
BES-246	Iris Straight		Needle, 0.5	10.8
BES-248	Iris Angled		0.5, 1.0	10.8

BES-246



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-250	Adson	1.0	12.0
BES-250-N	Adson	1.0	12.0
BES-252	Semkin	0.5, 1.0	14.0
BES-252-N	Semkin	0.5, 1.0	14.0
BES-254	Cushing	1.0, 1.5, 2.0	15.9
BES-256	Cushing	1.0, 1.5, 2.0	17.8



BES-258
Cushing Gursh




BES-260
Gerald Straight



BES-262
Gerald Angled



BES-264
Gerald Bayonet




Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Cushing	Tip Size (mm)	Length Size (cm)
BES-258	Gursh		1.0, 1.5, 2.0	17.8
BES-260	Gerald Straight		1.0, 1.5, 2.0	17.8
BES-262	Gerald Angled		1.0, 1.5, 2.0	17.8
BES-264	Gerald Bayonet		1.0, 1.5, 2.0	18.4

BES-
262
BES-
264



BES-266
Jansen Bayonet




BES-268
Jansen Bayonet



BES-270
Jansen Bayonet



BES-272
Malis



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm



Ref	Item Name	Jansen	Tip Size (mm)	Length Size (cm)
BES-	Bayonet		1.0, 1.5, 2.0	15.2
266	Jansen Bayonet		1.0, 1.5, 2.0	17.8
BES-	Jansen Bayonet		1.0, 1.5, 2.0	20.4
268	Malis		1.0, 1.5, 2.0	17.8

BES-

270

BES-

272



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-274	Scoville-Greenwood	1.0, 1.5	19.1
BES-276	Yasargil	1.0,	19.7
BES-278	Hardy Bayonet	1.5	21.6
BES-280	Hardy Bayonet Angled	1.0,	21.6
		1.5	
		1.0,	
		1.5	

BES-282
Jeweler # 4



BES-284
Jeweler # 7



BES-286
Adson



BES-288
Potts-Smith



BES-292
Potts-Smith




BES-294
Jansen Bayonet



BES-296
Jansen Bayonet



BES-302
Jansen Bayonet UP



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name Jeweller	Tip Size (mm)	Length Size (cm)
BES-	#4	Needle	11.5
282	Jeweller #7	Needle	11.5
BES-	Adson	1.0 2.0	12.0
284	Potts-Smith		17.8
BES-	Potts-Smith	2 .	21.0
286 0	Potts-Smith	0	24.1
BES-	Jansen Bayonet	2 .	13.3
288 2	Jansen Bayonet	0	15.9
BES-	Jansen Bayonet	2 .	19.1
298 8	Jansen Bayonet	0	21.0
BES-	Jansen Bayonet UP	2 .	19.1
290 0	Jansen Bayonet UP	0	21.0

BES- 2 .
3 0 2 0
BES- 2 .



BES-306-IR
Hardy Bayonet



BES-308-IR
Hardy Bayonet



BES-308-IR-N
Jansen Bayonet




BES-310-IR
Scoville-Greenwood



BES-312-IR
Gerald Bayonet



BES-314-IR
Straight



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-306-IR	Hardy Bayonet	0.5, 1.0, 1.5	21.6
BES-308-IR	Jansen Bayonet	0.5, 1.0, 1.5	19.1
BES-308-IR-N	Jansen Bayonet	0.5, 1.0, 1.5	19.1
BES-310-IR	Scoville-Greenwood	0.5, 1.0, 1.5	19.7
BES-312-IR	Gerald Bayonet	0.5, 1.0, 1.5	18.4
BES-314-IR	Straight	0.5, 1.0, 1.5	15.9

Note: More Irrigation Bipolar Forceps are available use IR extension to any regular bipolar forceps.



BES-316-NS-N
McPherson



BES-316-NS
McPherson



BES-318-NS
McPherson Angle



BES-320-NS
Jeweler # 5




BES-322-NS
Adson



BES-324-NS
Adson



BES-326-NS
Semkin



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP-** Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-316-NS	McPherson	0.5	10.2
BES-316-NS-N	McPherson	0.5	10.2
BES-318-NS	McPherson Angled	0.5	10.2
BES-320-NS	Jeweller #5 BES-322-	Needle, 0.5	11.5
NS Jeweller #4		Needle, 0.5	11.5
BES-324-NS	Adson	1.0	12.0
BES-326-NS	Semkin	0.5, 1.0	14.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.



BES-328-NS
Cushing



BES-332-NS
Gerald Bayonet



BES-334-NS
Jansen Bayonet




BES-336-NS
Scoville-Greenwood



BES-338-NS
Yasargil



BES-340-NS
Hardy Bayonet



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP-** Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-328-NS	Cushing	1.0, 1.5, 2.0	15.9
NS Cushing	BES-330-NS	1.0, 1.5, 2.0	17.8
Gerald Bayonet	BES-332-NS	1.0, 1.5	18.4
Jansen Bayonet	BES-334-NS	1.0, 1.5	15.9
BES-336-NS	Scoville-Greenwood	1.0, 1.5	19.1
BES-338-NS	Yasargil	1.0, 1.5	19.7
BES-340-NS	Hardy Bayonet	1.0, 1.5	21.6

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.

BES-341-HS
Jeweler # 4




BES-342-HS
Semkin




BES-343-HS
Jansen Bayonet



BES-344-HS
Scoville-Greenwood



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-341-HS	Jeweller #4	0.5,	11.5
HS	Semkin	0.	14.0
BES-343-HS	Jansen Bayonet	7	19.1
344-HS	Scoville-Greenwood	0.	19.7

7
1.
5

BES-345-HS
Jansen Bayonet



BES-346-HS
Jansen Bayonet




BES-347-HS
Potts-Smith



BES-348-HS
Potts-Smith



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP**- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-345-HS	Jansen Bayonet	2.0	19.0
BES-346-HS	Jansen Bayonet	2.	21.0
BES-347-HS	Potts-Smith BES-	0	21.0
348-HS	Potts-Smith	2.	24.0

0

2.

0



Tuner Warwick Scissors (STR)
BES-350



Tuner Warwick Scissors (CRV)
BES-352



Dressing Scissors (STR)
BES-358



Dressing Scissors (CRV)
BES-360



Wilson Hey Forceps (STR)
BES-366



Wilson Hey Forceps (CRV)
BES-368



Artery Forceps (STR)
BES-374



Artery Forceps (CRV)
BES-376



Abbey Needle
BES-382



Richie Forceps
BES-384

ENGLISH MONOPOLAR MALE 4.8MM CONNECTOR

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-350	Tuner Warwick Scissors (STR)	-	17.8
BES-352	Tuner Warwick Scissors (CRV)	-	17.8
BES-358	Tuner Warwick Scissors (STR)	-	23.0
BES-360	Tuner Warwick Scissors (CRV)	-	23.0
BES-366	Dressing Scissors (STR)	-	17.8
BES-368	Dressing Scissors (CRV)	-	17.8
BES-374	(CRV) Dressing Scissors	-	23.0
BES-376	(STR) Dressing Scissors	-	23.0
BES-366	Wilson Hey Forceps (STR)	-	17.8
BES-368	Wilson Hey Forceps (CRV)	-	17.8
BES-374	(CRV) Wilson Hey Forceps	-	23.0
BES-376	(STR) Wilson Hey Forceps	-	23.0
BES-374	Artery Forceps (STR)	-	17.8
BES-376	Artery Forceps (CRV)	-	17.8
BES-374	Artery Forceps (STR)	-	23.0
BES-376	Artery Forceps (CRV)	-	23.0
BES-382	Abbey Needle	-	17.5
BES-384	Richie Forceps	-	17.8
BES-384	Richie Forceps	-	20.0
BES-384	Richie Forceps	-	23.0

EUROPEAN MONOPOLAR FEMALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-390	Tuner Warwick Scissors (STR)	-	17.8
BES-392	Tuner Warwick Scissors (CRV)	-	17.8
BES-392	Tuner Warwick Scissors (STR)	-	23.0
BES-392	Tuner Warwick Scissors (CRV)	-	23.0
BES-394	Dressing Scissors (STR)	-	17.8
BES-396	Dressing Scissors (CRV)	-	17.8
BES-396	(CRV) Dressing Scissors	-	23.0
BES-400	(STR) Dressing Scissors	-	23.0
BES-400	Wilson Hey Forceps (STR)	-	17.8
BES-402	Wilson Hey Forceps (CRV)	-	17.8
BES-402	(CRV) Wilson Hey Forceps	-	23.0
BES-402	(STR) Wilson Hey Forceps	-	23.0
BES-410	Artery Forceps (STR)	-	17.8
BES-412	Artery Forceps (CRV)	-	17.8
BES-412	Artery Forceps (STR)	-	23.0
BES-412	Artery Forceps (CRV)	-	23.0
BES-420	Abbey Needle	-	17.5
BES-422	Richie Forceps	-	17.8
BES-422	Richie Forceps	-	20.0
BES-422	Richie Forceps	-	23.0



Adson
BES-430



McINDOE
BES-432



Gillies
BES-434



Dissecting
BES-442



Gerald
BES-450



Adson Bayonet
BES-452



Adson Bayonet UP
BES-456



Remington Hobb
BES-462

ENGLISH MONOPOLAR MALE 4.8MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
BES-430	Adson BES-	1.0 1.0,	12.0
432	McINDOE BES-	2.0 -	15.0
434	Gillies BES-436		15.0
	Gillies		17.8
BES-438	Gillies BES-	-	20.0
440	Gillies BES-442	1.0,	24.0
	Dissecting BES-444	2.0	15.0
	Dissecting	1.0,	17.8
BES-446	Dissecting BES-	2.0,	21.0
448	Dissecting BES-450	2.0	24.1
	Gerald 1x2 BES-452	1.0,	17.8
	Adson Bayonet	2.0 -	16.0
BES-454	Adson Bayonet	1.0, 2.0	19.1
BES-456	Adson Bayonet Angled UP	1.0, 2.0	16.0
BES-458	Adson Bayonet Angled UP	1.0, 2.0	19.1
BES-460	Adson Bayonet Angled UP	1.0, 2.0	21.0
BES-462	Remington Hobb	1.0, 2.0	24.1

EUROPEAN MONOPOLAR FEMALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
BES-364	Adson BES-	1.0 1.0,	12.0
366	McINDOE BES-	2.0 -	15.0
368	Gillies BES-370		15.0
	Gillies		17.8
BES-372	Gillies BES-	-	20.0
474	Gillies BES-476	1.0,	24.0
	Dissecting BES-478	2.0	15.0
	Dissecting	1.0,	17.8
BES-480	Dissecting BES-	2.0,	21.0
482	Dissecting BES-484	2.0	24.1
	Gerald 1x2 BES-486	1.0,	17.8
	Adson Bayonet	2.0 -	16.0
BES-488	Adson Bayonet	1.0, 2.0	19.1
BES-490	Adson Bayonet Angled UP	1.0, 2.0	16.0
BES-492	Adson Bayonet Angled UP	1.0, 2.0	19.1
BES-494	Adson Bayonet Angled UP	1.0, 2.0	21.0
BES-496	Remington Hobb	1.0, 2.0	24.1

ENGLISH MONOPOLAR MALE 4.0MM CONNECTOR


Ref	Item Name	Tip Size	Length Size
BES-498	Adson BES-	1.0 1.0,	12.0
500	McINDOE BES-	2.0 -	15.0
502	Gillies BES-504		15.0
	Gillies		17.8
BES-506	Gillies BES-	-	20.0
508	Gillies BES-510	1.0,	24.0
	Dissecting BES-512	2.0	15.0
	Dissecting	1.0,	17.8
BES-514	Dissecting BES-	2.0,	21.0
516	Dissecting BES-518	2.0	24.1
	Gerald 1x2 BES-520	1.0,	17.8
	Adson Bayonet	2.0 -	16.0
BES-522	Adson Bayonet	1.0, 2.0	19.1
BES-524	Adson Bayonet Angled UP	1.0, 2.0	16.0
BES-526	Adson Bayonet Angled UP	1.0, 2.0	19.1
BES-528	Adson Bayonet Angled UP	1.0, 2.0	21.0
BES-530	Remington Hobb	1.0, 2.0	24.1

MONOPOLAR WITHOUT CONNECTOR

Ref	Item Name	Tip Size	Length Size
BES-532	Adson BES-	1.0 1.0,	12.0
534	McINDOE BES-	2.0 -	15.0
536	Gillies BES-538		15.0
	Gillies		17.8
BES-540	Gillies BES-	-	20.0
542	Gillies BES-544	1.0,	24.0
	Dissecting BES-546	2.0	15.0
	Dissecting	1.0,	17.8
BES-548	Dissecting BES-	2.0,	21.0
550	Dissecting BES-552	2.0	24.1
	Gerald 1x2 BES-554	1.0,	17.8
	Adson Bayonet	2.0 -	16.0
BES-556	Adson Bayonet	1.0, 2.0	19.1
BES-558	Adson Bayonet Angled UP	1.0, 2.0	16.0
BES-560	Adson Bayonet Angled UP	1.0, 2.0	19.1
BES-562	Adson Bayonet Angled UP	1.0, 2.0	21.0
BES-564	Remington Hobb	1.0, 2.0	24.1



BES-570-DP
Jeweler # 5



BES-572-DP
Straight



BES-574-DP
Bayonet




BES-576-DP
Straight



BES-578-DP
Straight Angled



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: **TIP-** Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
BES-570-DP	Jeweller #5 With 3m Cable	0.5	11.5
BES-572-DP	Straight With 3m Cable	1.0,	18.0
BES-574-DP	Bayonet With 3m Cable	1.5	18.0
BES-576-DP	Straight	1.0,	18.0
BES-578-DP	Bayonet	1.0, 1.5	20.0

Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.
1.0
1.5



Tip size can be changed as per your requirement, just add the following extensions
fit tip size at the end of Ref #: TIP- Needle, 0.5, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	length Size (cm)
BES-580-DP	US. Gerald	1.0, 1.5	17.8
BES-582-DP	US. Jansen Bayonet Irrigation	1.0,	19.1
BES-584-DP	US. Jansen Bayonet BES-586-	1.5	17.8
DP US. Jeweller #5		1.0,	11.5
BES-588-DP	US. Adson	1.5, 0.5	12.1

Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.



Graves Vaginal Speculum
With Smoke Evacuation
Tube Small BES-600
Medium BES-602
Large BES-604

Disposable Tube



Small BES-606
Medium BES-608
Large BES-610



Pederson Vaginal Speculum
With Smoke Evacuation
Tube Small BES-612
Medium BES-614
Large BES-616

Disposable Tube



Small BES-618
Medium BES-620
Large BES-622



VU-Max Grave Vaginal Speculum
With Smoke Evacuation Tube
Medium **BES-624**
Large **BES-626**

Disposable Tube



Medium **BES-628**
Large **BES-630**



VU-Max Pederson Vaginal Speculum
With Smoke Evacuation Tube
Medium **BES-632**
Large **BES-634**

Disposable Tube



Medium **BES-636**
Large **BES-638**



Collin Vaginal Speculum
With Smoke Evacuation
Tube Small BES-640
Medium BES-642
Large BES-644



CUSCO Vaginal Speculum
With Smoke Evacuation
Tube Small BES-646
Medium BES-648
Large BES-650



VU-More Graves Vaginal
Speculum *With Smoke Evacuation
Tube* Medium BES-652
Large BES-654

Disposable Tube



Medium BES-656
Large BES-658



4-Blades Speculum
*With Smoke Evacuation
Tube* Medium BES-660
Large BES-662



Weisman Vaginal Speculum
Right Side Open
With Smoke Evacuation Tube
Medium BES-664
Large BES-666

Disposable Tube



Medium BES-668
Large BES-670



Weisman Vaginal Speculum
Left Side Open
With Smoke Evacuation Tube
Medium BES-672
Large BES-674

Disposable Tube



Medium BES-676
Large BES-678



Lateral Retractor
Wide Open BES-
680



Rigby Retractor
Wide Open
BES-682



Gelpi Retractor
Wide Open
BES-684



Kogan Speculum
3mm BES-686
4mm BES-688
5mm BES-690



Lateral Retractor
21cm 22mmx75mm
BES-692



Schroeder
Tanaculum 24cm
BES-694



Schroeder
Tanaculum 24cm
BES-696



Single Tooth
Tanaculum 24cm
BES-698



Campion
Forceps 2cm
BES-700



Spong
Forceps 2cm
BES-702

Dressing Forceps BES-704 21cm
BES-706 25cm



Tissue Forceps BES-708 21cm
BES-710 25cm



Iris Hook BES-712 25cm



Emmet Hook BES-714 25cm



Emmet Hook Straight BES-716 25cm



Emmet Hook Two Prong BES-718 25cm



Emmet Hook Three Prong BES-720 25cm





Graves Vaginal Speculum
With Smoke Evacuation
Tube Small BES-750
Medium BES-752
Large BES-754



Pederson Vaginal Speculum
With Smoke Evacuation
Tube Small BES-756
Medium BES-758
Large BES-760



Collin Vaginal Speculum
With Smoke Evacuation
Tube Small BES-762
Medium BES-764
Large BES-766



Lateral Vaginal Speculum
BES-766



Kogan Endocervical Speculum



3mm BES-770



4mm BES-772



5mm BES-774



Rigby Retractor
BES-778



Gelpi Retractor
BES-776



Schroeder Tenaculum
BES-780



Iris Hook
BES-782



Emmet Hook
BES-784



Emmet Hook Straight
BES-786



Emmet Hook Two Prong
BES-788



Emmet Hook Three Prong
BES-790

Without Suction



L-Hook Electrode 36cm
BES-800



J-Hook Electrode 36cm
BES-802



Needle Electrode 36cm
BES-804



Spectula Electrode 36cm
BES-806



Ball Electrode 36cm
BES-808



With Suction



L-Hook Electrode 36cm
BES-810



J-Hook Electrode 36cm
BES-812



Needle Electrode 36cm
BES-814




Spectula Electrode 36cm
BES-816



Ball Electrode 36cm
BES-818



Dagger 12xm
BES-820



Dagger 18cm
BES-822



Hook 12cm
BES-824



Hook 18cm
BES-826



Hockey 12cm
BES-828



Hockey 18cm
BES-830



Frazier Suction Tube 11FR 13.5cm
BES-832



Frazier Suction Tube 11FR 20.5cm
BES-834



Frazier Suction Tube 11FR 40.6cm
BES-836



Ball Electrode 36cm
BES-818



Rocker Switch Pencil (Single Use)
3m Cable BES-854



Foot Control Pencil (Single Use) 3m Cable BES-856
5m Cable BES-858



Finger Switch Pencil (Reusable)
3m Silicon Cable BES-860



Electrode Handle (Reusable)
3m Silicon Cable BES-862



Finger Switch Pencil (Single Use) 3m Cable BES-850
5m Cable BES-852



Blade Electrode

Disposable: 381-900-S Reusable: 381-900-R



Blade Electrode Crv

Disposable: 381-902-S Reusable: 381-902-R



Needle Electrode

Disposable: 381-904-S Reusable: 381-904-R



Needle Electrode Crv

Disposable: 381-906-S Reusable: 381-906-R



Ball Electrode 2mm

Disposable: 381-908-S Reusable: 381-908-R



Ball Electrode Crv 2mm

Disposable: 381-910-S Reusable: 381-910-R



Ball Electrode 3mm

Disposable: 381-912-S Reusable: 381-912-R



Ball Electrode 3mm Crv

Disposable: 381-914-S Reusable: 381-914-R



Ball Electrode 4mm

Disposable: 381-916-S Reusable: 381-916-R



Ball Electrode 4mm Crv

Disposable: 381-918-S Reusable: 381-918-R



Ball Electrode 5mm

Disposable: 381-920-S Reusable: 381-920-R



Ball Electrode 5mm Crv

Disposable: 381-922-S Reusable: 381-922-R



Blade Electrode Shaft 13cm

Disposable: 381-924-S Reusable: 381-924-R



Needle Electrode Shaft 13cm

Disposable: 381-926-S Reusable: 381-926-R



Ball Electrode 3mm Shaft 13cm

Disposable: 381-928-S Reusable: 381-928-R



Ball Electrode 5mm Shaft 13cm

Disposable: 381-930-S Reusable: 381-930-R



Conization Electrode

Disposable: 381-932-S Reusable: 381-932-R



Circle Loop Electrode 4mm

Disposable: 381-934-S Reusable: 381-934-R



Circle Loop Electrode 6mm

Disposable: 381-936-S Reusable: 381-936-R



Circle Loop Electrode 8mm

Disposable: 381-938-S Reusable: 381-938-R



Diamond Electrode

Disposable: 381-940-S Reusable: 381-940-R



Fine Needle Electrode

Disposable: 381-942-S Reusable: 381-942-R



					
Size	1.0x0.4cm	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm
Disposable	381-944-S	381-946-S	381-948-S	381-950-S	381-952-S
Reusable	381-944-R	381-946-R	381-948-R	381-950-R	381-952-R



					
Size	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm	1.0x1.5cm
Disposable	381-954-S	381-956-S	381-958-S	381-960-S	381-962-S
Reusable	381-954-R	381-956-R	381-958-R	381-960-R	381-962-R



						
Size	1.5x0.5cm	1.5x0.7cm	1.5x0.8cm	1.5x1.0cm	1.5x1.2cm	1.5x1.5cm
Disposable	381-964-S	381-966-S	381-968-S	381-970-S	381-972-S	381-974-S
Reusable	381-964-R	381-966-R	381-968-R	381-970-R	381-972-R	381-974-R



						
Size	2.0x0.8cm	2.0x1.0cm	2.0x1.2cm	2.0x1.3cm	2.0x1.5cm	2.0x2.0cm
Disposable	381-976-S	381-978-S	381-980-S	381-982-S	381-984-S	381-986-S
Reusable	381-976-R	381-978-R	381-980-R	381-982-R	381-984-R	381-986-R



Blade Electrode

Disposable: 381-1000-S Reusable: 381-1000-R



Blade Electrode Shaft 13cm

Disposable: 381-1024-S Reusable: 381-1024-R



Blade Electrode Crv

Disposable: 381-1002-S Reusable: 381-1002-R



Needle Electrode Shaft 13cm

Disposable: 381-1026-S Reusable: 381-1026-R



Needle Electrode

Disposable: 381-1004-S Reusable: 381-1004-R



Ball Electrode Shaft 3mm 13cm

Disposable: 381-1028-S Reusable: 381-1028-R



Needle Electrode Crv

Disposable: 381-1006-S Reusable: 381-1006-R



Ball Electrode 4mm Shaft 13cm

Disposable: 381-1030-S Reusable: 381-1030-R



Ball Electrode 2mm

Disposable: 381-1008-S Reusable: 381-1008-R



Ball Electrode 5mm 13cm

Disposable: 381-1032-S Reusable: 381-1032-R




Ball Electrode Crv 2mm

Disposable: 381-1010-S Reusable: 381-1010-R



Ball Electrode 3mm

Disposable: 381-1012-S Reusable: 381-1012-R



Circle Loop Electrode 4mm 0.7cm

Disposable: 381-1034-S Reusable: 381-1034-R



Ball Electrode 3mm Crv

Disposable: 381-1014-S Reusable: 381-1014-R




Circle Loop Electrode 6mm 0.7cm

Disposable: 381-936-S Reusable: 381-936-R



Ball Electrode 4mm

Disposable: 381-1016-S Reusable: 381-1016-R



Circle Loop Electrode 8mm 0.7cm

Disposable: 381-1038-S Reusable: 381-1038-R



Ball Electrode 4mm Crv

Disposable: 381-1018-S Reusable: 381-1018-R



Circle Loop Electrode 10mm 0.7cm

Disposable: 381-940-S Reusable: 381-940-R



Ball Electrode 5mm

Disposable: 381-1020-S Reusable: 381-1020-R



Plate Electrode Small

Disposable: 381-1042-S Reusable: 381-1042-R



Ball Electrode 5mm Crv

Disposable: 381-1022-S Reusable: 381-1022-R



Plate Electrode Large

Disposable: 381-1044-S Reusable: 381-1044-R



					
Size	1.0x0.4cm	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm
Disposable	381-1046-S	381-1048-S	381-1050-S	381-1052-S	381-1054-S
Reusable	381-1046-R	381-1048-R	381-1050-R	381-1052-R	381-1054-R



					
Size	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm	1.0x1.5cm
Disposable	381-1056-S	381-1058-S	381-1060-S	381-1062-S	381-1064-S
Reusable	381-1056-R	381-1058-R	381-1060-R	381-1062-R	381-1064-R



						
Size	1.5x0.5cm	1.5x0.7cm	1.5x0.8cm	1.5x1.0cm	1.5x1.2cm	1.5x1.5cm
Disposable	381-1066-S	381-1068-S	381-1070-S	381-1072-S	381-1074-S	381-1076-S
Reusable	381-1066-R	381-1068-R	381-1070-R	381-1072-R	381-1074-R	381-1076-R



						
Size	2.0x0.8cm	2.0x1.0cm	2.0x1.2cm	2.0x1.3cm	2.0x1.5cm	2.0x2.0cm
Disposable	381-1078-S	381-1080-S	381-1082-S	381-1084-S	381-1086-S	381-1088-S
Reusable	381-1078-R	381-1080-R	381-1082-R	381-1084-R	381-1086-R	381-1088-R



Bipolar Cable

Instruments Connector: USA 2-Pin Plug

Generator Connector: 4mm Banana Plug

Reusable Silicon 3m Cable BES-1100

Single Use 3m Cable BES-1102



Bipolar Cable

Instruments Connector: European Flat Plug

Generator Connector: 4mm Banana Plug

Reusable Silicon 3m Cable BES-1104

Single Use 3m Cable BES-1106



Bipolar Cable

Instruments Connector: USA 2-Pin Plug

Generator Connector: Pin Distance 29mm

Reusable Silicon 3m Cable BES-1112

Single Use 3m Cable BES-1114



Bipolar Cable

Instruments Connector: European Flat Plug

Generator Connector: Pin Distance 29mm

Reusable Silicon 3m Cable BES-1108

Single Use 3m Cable BES-1110



Bipolar Cable

Instruments Connector: Bipolar Artery Sealer

Generator Connector: 4mm Banana Plug

Reusable Silicon 3m Cable BES-1116



European Bipolar Cable

Instruments Connector: European Flat Plug

Generator Connector: Martin Type

Reusable Silicon 3m Cable BES-1118



US Bipolar Cable

Instruments Connector: USA 2-Pin Plug

Generator Connector: Martin Type

Reusable Silicon 3m Cable BES-1120



European Bipolar Cable

Instruments Connector: European Flat Plug

Generator Connector: Erbe Type



Reusable Silicon 3m Cable BES-1122

US Bipolar Cable

Instruments Connector: USA 2-Pin Plug

Generator Connector: Erbe Type



Bipolar Cable

Instruments Connector: Ethicon

Generator Connector: 4mm Banana Plug

Reusable Silicon 3m Cable BES-1126



Monopolar Cable
Instruments Connector: Screw Type
Generator Connector: Standard Bovie 8mm Plug



Reusable Silicon 3m Cable BES-1128
Monopolar Cable
Instruments Connector: 4mm Female Fitting
Generator Connector: Standard Bovie 8mm Plug



Monopolar Endoscopic Cable
Instruments Connector: 4mm Female Fitting
Generator Connector: 4mm Banana Plug
Reusable Silicon 3m Cable BES-1132



Diathermy Cable
Instruments Connector: 4mm Male Fitting
Generator Connector: 4mm Banana Plug
Reusable Silicon 3m Cable BES-1134



Diathermy Cable
Instruments Connector: 4.8mm Female Fitting
Generator Connector: 4mm Banana Plug
Reusable Silicon 3m Cable BES-1136



Diathermy Cable
Instruments Connector: 4.8mm Female Fitting
Generator Connector: Stan. Bovie 8mm Plug
Reusable Silicon 3m Cable BES-1138

BES-1150
BES-1152
BES-1154

BES-1156
BES-1158
Artery Sealer 6.5"
Artery Sealer 7.5"
Artery Sealer 8.5"



BES-1160
BES-1162
BES-1164

BES-1166
BES-1168
Non-Stick Artery Sealer 6.5"
Non-Stick Artery Sealer 7.5"
Non-Stick Artery Sealer 8.5"





BEST ELECTRO
SURGICAL



- 📍 Adalat Garh, Tehsil & District
Sialkot 51310 Pakistan +92-
+92-329-6320000
- ☎ info@bestelectrosurgical.com
- 🌐 www.bestelectrosurgical.com